




the best MRI experience







Stress free The best MRI experience

Our key value has always been, since the very beginning, Patients First. MROpen Evo provides superior comfort, greatly reducing the problem of claustrophobia and patient anxiety that often occurs in closed "tunnel" systems, causing concern and distress in both patients and MRI operators.

The Patient can walk into the scanner, can sit, lie slightly backward, lie horizontally or even stand. With no barrier with the surrounding environment, patients can see around them at all times or enjoy watching TV while comfortably undergoing an MRI procedure.

All of this while improving the clinical outcome by being able to scan the patient in the position of maximum symptom.




Stand out

Make your MRI centre truly unique

Thanks to its unique design MROpen Evo offers something different to both patients and clinicians.

The unique vertical opening accommodates larger, claustrophobic, elderly and young patients, or those who are in pain, with no barrier between the patient and the environment.

To cover specific needs, differentiate your services and increase productivity MROpen Evo is the perfect choice for your centre as a stand-alone scanner or associated with a traditional MRI.




Don't guess. See

Upright and multi-position imaging

MROpen Evo changes the approach to the clinical investigation using positional symptom-driven imaging.

New clinical advantages come from MROpen Evo's Upright and multi-position imaging capabilities of the Spine and Joints in Weight Bearing mode as well as functional studies of the Spine and Joints in flex-extension.

Changing point of view makes the difference.






SUPINE FSE T2
4 mm thk - 5'10"




SUPINE SE T1
4 mm thk - 4'24"



SUPINE FIR PD
4 mm thk - 5'19"



SUPINE FSE T2
4 mm thk - 4'09"




SITTING FSE T2
4 mm thk - 4'53"



STANDING FSE T2
5 mm thk - 1'19"



Cervical Spine




Routine
Flexion/Extension
Soft tissue neck
MRA carotids


Patient Positions:
Supine
Angle 20
Angle 40
Sitting



SUPINE FSE T2
3 mm thk - 5'40"



SUPINE GE T2
4 mm thk - 4'02"



FLEX FSE T2
4,5 mm thk - 1'29"



EXT FSE T2
4,5 mm thk - 1'49"

Thoracic Spine Lumbar Spine




Routine
Flexion/Extension
Weight Bearing

Patient Positions:
Supine
Angle 20
Angle 40
Angle 60
Sitting
Standing



Brain


IAC
TMJ
Pituitary gland
Orbits
MRA



Patient Positions:
Supine
Angle 20
Angle 40
Sitting



SUPINE SE T1
5 mm thk - 4'27"



SUPINE FSE T2
5 mm thk - 4'10"



SUPINE FLAIR T2
5 mm thk - 4'36"



SUPINE MRA
1,2 mm thk - 6'27"




Knee



Patient Positions:
Supine
Angle 20
Angle 40
Angle 60


Sitting
Standing
Bent knee



SUPINE SE T1
3 mm thk - 5'13"




SUPINE FIR PD
3 mm thk - 6'36"



SUPINE FSE T2
3 mm thk - 5'40"



SUPINE FSE T2
4,5 mm thk - 5'21"




STANDING FSE T2
4,5 mm thk - 2'13"




STANDING SE T1 -
4 mm thk - 2'41"




SUPINE FIR PD
4 mm thk - 5'00"



SUPINE GE T2
4 mm thk - 3'29"




SITTING SE T1
3 mm thk - 2'13"




SUPINE SE T1
4 mm thk - 2'28"

Extremities




Patient Positions:
Supine
Angle 40
Sitting
Standing




Patient Positions:
Supine
Angle 40
Sitting
Standing




Hip



Patient Positions:
Supine
Standing
Standing with hip rotation




SUPINE FSE PD
5,5 mm thk - 5'16"




STANDING SE T1
4,5 mm thk - 2'44"




STANDING FSE T2
5 mm thk - 1'38"




STANDING FSE T2
5 mm thk - 1'38"

Receiving Coils




KNEE
Channels: 2
Clinical Applications: Knee, Ankle, Foot



BODY SPINE
Channels: 4
Clinical Applications: L-Spine, T-Spine, Pelvis, Bilateral Hips, Sacroiliac Joints



WRIST-HAND
Channels: 2
Clinical Applications: Hand, Wrist, Elbow




MP-FLAT
Channels: 1
Clinical Applications: Upright/Sitting Hip Pelvic Floor, Multipurpose



SPINE FLAT
Channels: 4
Clinical Applications: Upright/Sitting, L-Spine, T-Spine, Sacroiliac Joints



CERVICAL SPINE
Channels: 2
Clinical Applications: C-Spine



MP-LOOP
Channels: 1
Clinical Applications: Shoulder, Dynamic studies, Multipurpose




SHOULDER
Channels: 2
Clinical Applications: Shoulder




WHOLE SPINE
Channels: 8
Clinical Applications: Spine



FLEX SMALL/LARGE
Channels: 1
Clinical Applications: Upright/Sitting Pelvis, Bilateral Hips



BRAIN
Channels: 2
Clinical Applications: Head - Brain



HEAD NECK
Channels: 3
Clinical Applications: Brain, MRA

The Helium-free MRI



UNIQUE DESIGN

U-shaped design delivers unparalleled patient comfort with no barrier between the patient and the environment.

UPRIGHT IMAGING

New approach towards pathologies of the Spine and Joints, by allowing earlier detection and more accurate diagnosis.

MgB₂ HELIUM-FREE TECHNOLOGY

Innovative superconducting material (Magnesium Diboride – MgB₂).
Environmentally friendly.
Quench-free safe operation.
Lower power consumption.

HIGHEST PATIENT ACCEPTANCE

Claustrophobic, anxious, larger, elderly or disabled patients, or those who are in pain, are able to undergo the examination in comfort and without anxiety.

MULTI-POSITION IMAGING

Imaging can be performed with the patient positioned appropriately according to symptoms: standing, sitting, bending or lying down.

SPECIAL STUDIES & RESEARCH

Unmatched versatility enables special MRI studies and clinical research: from dynamic studies, to imaging professional athletes with extreme ranges of flexion of the spine and joints.



Testimonials

Patients

"I'm a combat veteran from US Air Force and I'm suffering from PTSD. A closed bore MRI was not an option for me. In the MROpen Evo everything is different because you can sit and see all around you"



Over the years my claustrophobia increased and it left me dreading necessary MRI scans. After experiencing the MROpen that all changed.



I had a problem with my lumbar spine. After the examination in the MROpen scanner, being in the flexion and extension position as well I got to see what my real problem had been this whole time.



A child suffering from Asperger syndrome, Tourette syndrome, and ADHD has had many terrible experiences with MRIs and finds relief with the MRI system at an Innovative MRI Partners facility.

Testimonials

Professionals

Penny Gowland
Professor of Physics, Faculty of Science
Sir Peter Mansfield Imaging Centre
Nottingham University



Open scanners give a whole new dimension to MRI, which can be used to study the effects of gravity on structure and function, and also dynamic processes in the human body. This will provide new clinical information most obviously in MSK imaging but also in lung, GI and neuroimaging, and will also open up new avenues in experimental medicine research. Furthermore, the convenience of the scanner offers the chance for faster patient scanning with less risk of claustrophobia.

Alan Breen DC
PhD Professor of MSK Healthcare
AECC University College
Bournemouth (UK)



MROpen can be used both to study the service uses of such scanners and to conduct spinal biomechanics research that uses positional data (e.g. measuring and comparing disc height, volume and area in different positions).



MgB₂

The true Helium-free MRI innovation



SMART

MROpen EVO makes the best use of the most recent MgB₂ high-temperature superconductor and helium-free cooling solutions in an open-sky magnet, enabling an innovative approach to imaging and therapy.



SAFE

MgB₂ superconductor and helium-free cooling combined in a very stable magnet solution allow full control to the user, including switching on-off, and simplified start-up and maintenance.



GREEN

Much reduced use of helium and of rare-earth elements compared to other magnet solutions and full reliance on electricity to operate all help to protect our environment.



SIMPLE

Eliminating the need for gas refills or helium venting lines makes magnet installation and running much simpler than with conventional solutions.



ADVANCED

We allow you to do what you would expect from a conventional MRI system and in the meantime we allow you to think differently and provide innovative clinical solutions that others cannot, always for better patient outcomes and cost savings throughout cycle of care.

Technical Specification

Field Strength	0.5 T
Type	Superconductive “Cryogen-free” based on the technology of MgB2 (Magnesium Diboride)
Lateral Gap	56 cm 22 ³ / ₆₄ in
Digital Spectrometer	I-box spectrometer
Receiving channels	8
Gradients max intensity	20 mT/m
Rise time	0.6 msec (0 – 20 mT/m)
Slew rate	33 mT/m/msec
Receiving coils technology	Multi array cross tuned
Max patient weight	200 kg / 440 lb
Patient position	Laying, sitting, standing
Graphic user interface	MR-GUI Pro

SITING REQUIREMENTS

Total System Weight	28,000 kg / 61729 lb
Dimensions	170x200x200 cm (WxDxH) / 66 ⁵⁹ / ₆₄ x 78 ⁴⁷ / ₆₄ x 78 ⁴⁷ / ₆₄ in
Electronic Cabinets	3 x 250 kg / 3 x 551 lb
RF Shielding	Required
Power requirement	400V 3N ~50Hz 480V 3N ~50/60 Hz
Power Consumption	35 KVA

PULSED SEQUENCES

Spin Echo
Gradient Echo(*)
Multi Echo, PD-T2
Gradient Echo Dual Echo (HiCon) (*)

Fast Spin Echo

Fast Spin Echo (FSE T2, FSE T1, FSE PD)(*)
Fast Inversion Recovery (FIR PD, FIR T2)(*)
Fast Spin Echo Fluid Attenuated Inversion Recovery (Fast FLAIR)(*)
Metal T1/T2
HASTE
Diffusion Weigthed Imaging (DWI - HASTE)

Gradient Echo in Steady State

Balanced Steady State in Free Precession (GBASS 3D, GBASS 2D)
Gradient Echo in Steady State (GE 3D)
Gradient Echo in Steady State RF Spoiled (GE T1 3D)
Time Reversed Gradient Field Echo (EMIT 3D)
Dual Echo in Steady State (STSS 3D)

MRA

Angio 2D
Angio 3D

Fat Suppression

Gradient Echo Fat Water Separation (FWS)
Short time Inversion Recovery (GE-STIR, SE-STIR)

(*) Flow Compensation available





evo MROpen

ASG
SUPERCONDUCTORS

the best MRI experience

www.mropenevo.com

